

Fixi and Fox YANO

Wheeled Forearm Walker for Children



Fixi and FoxYano:

- Comes standard with 3 handle options to accommodate all grip types
- Dual action brake
- Foldable Design
- Airless Tires
- Powder coated aluminium frame
- Approved up to 100 kg/220 lbs SWL


ref.no.186.99...

Accessories:


- One handed brake (specify required hand)
- ref.no.280.99.70



Standard-handle



Ellipse-handle



Spade-handle

Type:	Fixi YANO	FoxYANO
Depth:	26.57"	26.57"
Width:	23.03"	23.03"
Height of forearm support:	35.62"-39.17"	38.97"-42. 71"
Adjustments in height:	5x	5x
Max. body weight:	100 kg/220 lbs	100 kg/220 lbs
Ref.no:	186.99.10	186.99.20