

Fixi and Fox

Rollator Wheeled Walker for Children



Fixi and Fox:

- Powder coated aluminium frame
- Dual action brake
- Foldable design
- Airless tires
- Approved up to 100 kg/220 lbs SWL

ref.no.186.00 ...

Accessories:

- Seat (black)
ref.no.286.65.25
- Back pad
ref.no.280.10.48



Type:	Fixi	Fox
Depth:	26.77"	26.77"
Width:	23.22"	23.22"
Height of handle:	24.01"-27.16"	27.16"-28.14"
Seat height (optional):	18.11"	18.11"
Seat width (optional):	14.96"	14.96"
Adjustments in height:	5x	5x
Max. body weight:	100 kg/220 lbs	100 kg/220 lbs
Ref. no.:	186.00.10	186.00.20



Technical changes and changes in colour reserved.

protected
by
patent

CE