

Fun-Kids

Crutches with open forearm cuff




Fun-Kids:

- Colorful, suitable for children
- 3 adjustments in height from handle to cuff
- Ergonomic handle
- Reflector in handle
- 10 adjustments in height
- Compression fitting
- Approved up to 100 kg/220 lbs SWL

ref.no.122.95

Accessories:

- Handle upholstery
- ref.no.160.01.01



Type:	Fun-Kids
Height from handle to cuff:	6.29"-8.26"
Height of handle:	19.68"-28.54"
Diameter in inches:	0.78"/0.66"
Ref.no.	122.95