

Pluto

Posterior Reverse Walker for Children



Pluto:

- 2 available sizes
- Width and height adjustable hand grips
- Return motion break
- Steerable front wheels
- Airless tires

ref.no.187.00. ..

- The Pluto reverse walker provides independent movement for children
- It supports the user from behind and promotes upright posture when standing and walking
- Supports therapeutic exercises for natural movement and to increase balance
- The braking system in the back wheels prevents unintentional backwards rolling

Accessories:

- Front-wheels with steerlock
ref.no.287.10.42
- Seat for Pluto - size 1 : up to 60 kg/133 lbs
ref.no.287.20.10



Type:	Size 1	Size 2
Length:	31.49"	33.85"
Width:	25.59"	26.37"
Height of handle:	18.11"-31.88"	21.65"-35.03"
Max. body weight:	60 kg/132 lbs	80 kg/176 lbs
Ref. no.:	187.00.10	187.00.20

Technical changes and changes in colour reserved.

