

Safe-In-Kids

Cruches with Forearm Cuff



Safe-In-Kids:

- Forearm cuff with hinge to increase mobility
- 3 adjustments in height from handle to cuff
- Ergonomic handle
- Reflector in handle
- 10 adjustments in height
- Compression fitting
- Solid ferrule with a high diameter
- Approved up to 100 kg/220 lbs SWL

ref.no.122.96.10



Type:	Safe-In-Kids
Height from handle to cuff:	6.29"-8.26"
Height of handle:	19.68"-28.54"
Diameter in inches:	0.78"/0.66"
Ref.no.	122.96.10