

# Walk-on

## Walking Frame with 2 Wheels for Children



**REBOTEC®**



### Walk-on:

- Powder coated aluminium frame
- Non Slipping Handles for safer use
- Foldable, height adjustable
- Approved up to 120 kg/265 lbs SWL

ref.no.175. ...

| Type:                  | XS             | S              | M              |
|------------------------|----------------|----------------|----------------|
| Depth:                 | 16.14"         | 16.14"         | 18.11"         |
| Width:                 | 22.24"         | 22.24"         | 24.60"         |
| Height of handle:      | 24.40"-27.75"  | 27.55"-31.69"  | 26.77"-34.44"  |
| Adjustments in height: | 5x             | 5x             | 10x            |
| Max. body weight:      | 120 kg/265 lbs | 120 kg/265 lbs | 120 kg/265 lbs |
| Color:                 | Multicolor     | Multicolor     | Charcoal       |
| Ref.no.                | 175.10         | 175.20         | 175.30         |

Technical changes and changes in colour reserved.

