

YANO Walker

Wheeled Forearm Walker



YANO Walker:

- Comes standard with 3 handle options to accommodate all grip types
- Height adjustable
- Powder coated frame
- 4 locking castors
- Approved up to 150 kg/330 lbs SWL


ref.no.188.99.10



Type:	YANO Walker
Depth:	28.34"
Width:	26.77"
Height of forearm support:	26.77"-45.66"
Adjustments in height:	10x
Max. body weight:	150 kg/330 lbs
Ref. no.:	188.99.10



Standard-handle



Ellipse-handle



Spade-handle